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OPERATION MANUAL

CONTENTS

1	DEVICE DESCRIPTION AND OPERATION.....	23
1.1	Function and Application Range	23
1.2	Technical Characteristics	24
1.3	Measuring Method	25
1.4	The Device Design General Data.....	26
2	MEASURING PROCEDURE.....	30
2.1	Operational Restrictions and Safety Measures	30
2.2	The Device Preparation for Operation	30
2.3	The Device Operation	31
2.3.1	Gamma Radiation Dose Power Measuring	31
2.3.2	Beta Particle Flux Density Measuring.....	31
2.3.3	Alpha Particle Flux Density Indication	33
2.3.4	Research and Checking of Items or Samples Contami- nated by Radioactive nuclides	34
3	TECHNICAL MAINTENANCE	35
3.1	Operation and Transportation Safety Measures	35
3.2	Technical Maintenance Procedure	35
4	TECHNICAL CERTIFICATE DATA	36
4.1	Package Contents	36
4.2	Service Life and Guarantee Liabilities	36
4.3	To be filled in by the manufacturing plant	37

1 DEVICE DESCRIPTION AND FUNCTIONING

1.1 Function and Application Range

Expert-M Dosimeter (hereinafter referred to as the device) is meant for photonic (roentgen and gamma) radiation (hereinafter referred to as gamma radiation dose power) ambient equivalent power measuring, for beta particle flux density measuring and for alpha particle flux density estimation, such this flux being emitted from contaminated surface. The device gives the opportunity to estimate beta, gamma and alpha radiating nuclides contamination of water, ground, foodstuff, agricultural crop and livestock farming product etc. samples and to perform operative search of radioactive emission sources or of contaminated objects.

The device maintains estimation of workplace, dwelling and ground radiation safety, radioactive contamination estimation as regards real objects, items, materials and samples including paper money and its packing. The device is usable; there are only two controls. Data is displayed on a four-digit indicator panel. Continuous measuring mode is used in the device and data showing measured values on the indicator panel is updated every 2 seconds; it is very handy for operative monitoring. The device has audio indication maintaining search of radiation sources via audio signaling frequency change.

The device power supply is maintained by a PP3 battery type battery (9 V).

1.2 Technical Characteristics

1.2.1 The device has the following technical characteristics:

Photonic radiation dose power measuring range in mcSv/hour (mcR/Hour)	0.10-100 (10-10000)
Photon power range in dose power level measuring, MeV	0.04-3.0
Energy dependence in measuring of dose power, %, within the limits	±30
Beta particle flux density measuring range, this flux being emitted from contaminated surface (Strontium-90 – Yttrium-90 particles), part./(cm ² .min)	3-10000
Lower limit of registered beta radiation energy, not exceeding, in MeV	0.05
Basic error, in %, within the limits	±25
Average value of the natural radiation background in γ-mode, mcSv/hour (mcR/Hour)	0.08 (8)
Indication range for alpha radiation flux density, this flux being emitted from contaminated surface (Plutonium-239 particles), part./(cm ² .min)	10 ² -10 ⁵ (или 10 ³ -10 ⁶)
Lower limit of registered alpha radiation energy, not exceeding, in MeV	3.0
Anisotropy of sensitivity for Cesium-137 and Americium-241 in spatial angle 4π, in %, within the limits	±40
Operating mode timing, in minutes, not exceeding	1
Shift/ data timing, from	2/40
Continuous function duration (for measuring performed at the natural radiation background level), in hours, not less than: - a PP3 battery (9 V)	50
Operational conditions:	

- Temperature, °C - Humidity, %	From -20 up to +40 Up to 90 at T 30°C
Dimensions, in mm	240x78x75
Weight (including accumulator), in g, not exceeding	500(550)

1.3 Measuring Method

An end-located gas-discharge counter with a thin detector window is applied in the device for radiation detection. Gamma quantum, beta or alpha particle flux is transformed by the detector into electric signal sequence. These signals are formed as per their duration and amplitude and then they are transmitted to registration and display circuit.


On the monitor average frequency of impulses in dose power units is displayed permanently in "γ" mode or beta and alpha particle flux density is displayed continuously in "β" mode. If value of irradiation level is constant, then data settling time is about 40 seconds. Indication on the display is updated every 2 seconds; so in case positional relationship between the source and the device has constant geometry, it is necessary to monitor the indications during 10-20 seconds after the data settling and to determine their average value in order to improve measuring accuracy. If irradiation level or distance between the source and the device change, indication starts changing immediately and the new value of indication settles after about 40 seconds. Change of indication every 2 seconds reduces the irradiation source search time.

The device generates a sound signal in response to each counter impulse; this function gives the opportunity to detect the radiation level change tendency evaluating average frequency of sound signals by hearing.

1.4 The Device Design General Data

1.4.1 The following details are placed on the device top panel (Picture 1):

- Digital display;
- Controls – two switches.

Picture 1

1.4.2 The display indicates the measured parameter in the digital mode (Pictures 1a, 1b):


Mode γ

Result – 0,56 mkSv/h

Picture 1a

Mode β

Result – 256 particles/cm² min

Picture 1b


If the power supply source is discharged three points are indicated on the display continuously (Picture 2).

Power supply source discharge

Picture 2

If the measured value range upper limit is overcome, "1" is indicated in the senior digit on the display; all the other digits at that are not indicated (Picture 3).

Upper limit excess

Picture 3

1.4.3 Controls on the top cover (Picture 1):

- Power supply switch is placed on the left and it has two positions: "ON" and "OFF";
- Mode switch is placed on the right and it has two positions:
 - " γ " – gamma radiation dose power measuring mode;
 - " β " – beta particle flux density measuring mode or alpha particle flux density indication mode – related to the measuring method applied.

1.4.4 The following details are placed on the back panel of the device (Picture 4):


- Detachable detector screen;
- Detachable cover of power supply compartment.

1.4.5 Detachable screen is meant for energy dependence data equalizing in the course of various energy gamma radiation registration as well as for the detector screening against beta and alpha particles in " γ " measuring mode.

The screen is fixed by a latch. It is necessary to push the protruding latch up at the point indicated in Picture 4 in order to detach

the screen. When fixing the screen it is required to push the same latch down.

1.4.6 The detachable cover of the power supply compartment maintains access to the power supply battery. In order to detach the cover it is necessary to push it up at the points indicated in Picture 4; in order to put the cover back it is required to push it down.

Picture 4

2 MEASURING PROCEDURE

2.1 Operational Restrictions and Safety Measures

2.1.1 Protect the detector inlet window covered by a net because its damage leads to the counter breakage.

2.1.2 It is prohibited to open the sealed compartment of the device in order to prevent contacting with high voltage in the detector power supply and breakage of the indicator circuit components.

2.1.3 Keep the power supply compartment and power supply contacts clean.

2.1.4 Replace discharge power supply batteries in good time.

2.1.5 If radioactive substance is found on the device case, its background indications may increase. Check it by measuring the device background indications at another spot or in another room.

2.2 The Device Preparation for Operation

It is necessary to perform the following actions in order to prepare the device for operation:

- Power supply switch shall be put into OFF position;
- The power supply compartment cover shall be removed; see Picture 4 where the spots to be pushed up are shown;
- a PP3 battery shall be put in observing polarity of connection;
- The power supply compartment cover shall be closed. The device is ready for operation 1 minute after switching it on by putting the power supply switch into ON position;

2.3 The Device Operation

2.3.1 Gamma Radiation Dose Power Measuring

It is required to implement the following actions for gamma radiation dose power measuring indoors or outdoors:

- The working surface of the detector shall be covered by the screen pursuant to Para.4.1;
- The device shall be switched on;
- "γ" operation mode shall be set;
- The device shall be placed at the distance of at least 1 m from the floor (ground) and from any surrounding objects;
- N dose power value in mSv/ hour shall be determined after 40 seconds (or more).

If more accurate measuring is required, n of N_{γ_i} device indication counts shall be registered and arithmetic average of N_{γ} dose power values in mSv/ hour shall be calculated using the following formula:

$$N_{\gamma} = \sum_{i=1}^n N_{\gamma_i} / n \quad (1)$$


Recommended number of measuring is from 5 up to 10. If it is required to measure small values of the dose power more accurately (within the range from 0.1 up to 1.0 mSv/ hour), it is necessary to deduct the value of the counter intrinsic radiation background correction from the device readings. This correction is equal to 0.5 of the device intrinsic radiation background value and at the average it makes 0.04 mSv/ hour.

2.3.2 Beta Particle Flux Density Measuring

In order to measure density of beta particle flux emitted from the surface of items contaminated by beta emitting radionuclides the following actions shall be performed:

- The detector working surface shall be covered by a screen;
- "β" operational mode shall be set;
- The device shall be switched on;

- The device shall be placed at the distance of 3-5 mm from the measured object surface (see Picture 5); 40 seconds after that n of counts completed by the device with covered $N_{3\gamma I}$ detector shall be registered (n from 5 up to 10) in part./ (cm²·min.), which is conditioned by background gamma radiation and by gamma radiation of the item being checked;

Picture 5

- Arithmetic average value of the device background readings showing its gamma radiation in "β" $N_{3\gamma}$ mode shall be calculated in part./ (cm²·min.) pursuant to the following formula:

$$N_{3\gamma} = \sum_{i=1}^n N_{3\gamma_i} / n \quad (2)$$

- The detector screen shall be detached and the device shall be placed on the same checking spot at the distance of 3-4 mm from the surface of the measured object again;

- 40 seconds later m of the device indication counts (m from 5 up to 10) with the open N_{oi} detector shall be performed and arithmetic average value of the N_o device summarized readings based on beta and gamma radiation shall be calculated in part./($\text{cm}^2 \cdot \text{min.}$), pursuant to the following formula:

$$N_o = \sum_{i=1}^m N_{oi} / m \quad (3)$$

N_β beta particle flux density shall be determined in part./($\text{cm}^2 \cdot \text{min.}$), this flux being emitted by the measured object contaminated by beta radiating nuclides pursuant to the following formula:

$$N_\beta = N_o - N_{\gamma} \quad (4)$$

2.3.3 Alpha Particle Flux Density Indication

The following actions are required to estimate (indicate) alpha radiation flux density, this flux being emitted from the surface of objects contaminated by alpha emitting nuclides:

- The detector screen shall be detached;
- "β" operational mode shall be set;
- The device shall be switched on;
- The device shall be placed at the distance of 3-5 mm from the measured object surface (see Picture 5), and 40 seconds after that n of summarized N_{ci} device indication readings based on alpha, beta and gamma radiation (n from 5 up to 10) shall be registered in part./($\text{cm}^2 \cdot \text{min.}$);

- Average arithmetic value of those summarized N_c readings shall be determined in part./($\text{cm}^2 \cdot \text{min.}$), pursuant to the following formula:

$$N_c = \sum_{i=1}^n N_{ci} / n \quad (5)$$

- The device shall be located at the distance of 30-50 mm from the surface of the object being studied (in this way maintaining ab-

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- Average arithmetic value of those summarized N_c readings shall be determined in part./($\text{cm}^2 \cdot \text{min.}$), pursuant to the following formula:

$$N_{cl} = \sum_{i=1}^m N_{cli} / m \quad (6)$$

N_{α} alpha particle flux density in part./ $\text{cm}^2 \cdot \text{min.}$, this flux being emitted from the studied object surface, shall be determined pursuant to the following formula:

$$N_{\alpha} = (N_c - N_{cl}) \cdot K_{\infty} \quad (7)$$

where K_{∞} is coefficient of sensitivity to alpha radiation (approximately 10^3).

2.3.4 Research and Checking Of Items or Samples Contaminated By Radioactive Nuclides

Research and Checking of items or samples to detect radioactive contamination shall be performed in order to detect certain objects (for example, construction materials, paper money etc.) or samples (of ground, foodstuff, agricultural products etc.) contaminated by radionuclides. These works shall result in sorting of checked objects or kinds of products in accordance with the permissible normative levels of radioactive contamination accepted for the latter.

Measurements related to the said works shall take into account specificity and physical characteristics of the checked objects as well as the tasks appearing while arranging such checking.

3 TECHNICAL MAINTENANCE

3.1 Operation and Transportation Safety Measures

3.1.1 Before starting the device operation it is necessary to read this Operation Manual attentively. It is prohibited to open or to repair the device for the Customer operating it because there is high voltage maintaining the counter power supply inside it, namely, about 400 V. So if repair is required it is expedient to send the device to the manufacturer or to employ persons who know the device circuit and design well.

But in order to maintain radiation hygiene it is recommended to protect the device against destruction or loss and to open the holder cover only in the course of maintenance works.

3.1.2 While preparing the device for transportation by any postal service or by means of transport **it is required to switch it off and to take the autonomous power supply source (the battery) out of it in order to avoid accidental triggering of the device which generates sound signals at that.**

3.2 Technical Maintenance Procedure

The device technical maintenance is performed in order to maintain its operating capacity during the service life term; it shall be performed by persons working with the device taking into account safety measures in accordance with Para.3.1.

Preventive works performed in the course of technical maintenance include package contents checking, visual examination of the device appearance and its serviceability checking. The device package content checking is to be completed by comparing its compliance with Para. 4.1. While examining the device appearance it is required to check absence of chippings and cracks on the device case, of mica turbidity on the counter inlet window as well as safe securing of switches.

4 TECHNICAL CERTIFICATE DATA

4.1 Package Contents

No.	Description	Number of pieces
1	EXPERT-M Dosimeter	1
2	PP3 battery	1
3	Operation Manual	1

4.2 Service Life and Guarantee Liabilities

Average device service life until complete overhaul is at least 6 years. Upon expiration of the said term further utilization of the device is possible after complete overhaul performed by the Manufacturer. The set service life after complete overhaul is fixed by the Manufacturer in the warranty service coupon issued for repair.

The Manufacturing Plant address is given in Para. 4.3.

The Manufacturing Plant guarantees serviceability of the device provided the Customer follows operation rules described herein.

The device guarantee service life is 24 months from the date of purchase (if it was sold to the Customer through the distribution network). The Manufacturing Plant performs guarantee, post-guarantee and complete overhaul. The Manufacturing Plant address is given in Para. 4.3.

The device guarantee repair term is not included into the set guarantee service life term. Claims are not accepted and guarantee repair is not performed in the following cases: the Customer's negligent device handling, the detector inlet window, indicator and the device case damage or if the device seals are broken.

4.3 To Be Filled In By the Manufacturing Plant

EXPERT-M Dosimeter, Manufacturer's Number _____ ,
was acknowledged serviceable.

Manufacturing Date _____

Seal

The Manufacturing Plant Address:

LLC "Axelbant"

125367, Moscow, Gabrichevskogo str., 10-4-396.

To be filled in by the trading Company:

Manufacturing Date _____

Seller _____

Stamp